

The Pathology Center
8303 Dodge Street
Omaha, NE 68114
402-354-4541
888-432-8980
Fax 402-354-8806

October 31, 2008

Update on Vitamin D, 25-hydroxy (25-OH Vitamin D)

Due to increased requests for Vitamin D 25-hydroxy, specimens are now run Monday, Wednesday, and Friday. Specimens are stable up to 5 days refrigerated and do not need to be protected from light.

Update on Stool Giardia and Cryptosporidium Tests

A new method for giardia antigen testing and cryptosporidium antigen testing has been implemented. Specimens for these 2 tests must be placed in a formalin vial (pink top) within 30 minutes of collection. Once the specimen is placed in formalin, it is stable for 2 weeks room temperature or refrigerated.

Please contact The Pathology Center at 402-354-4541 or 888-432- 8980 if there are any questions.