

8303 Dodge Street Omaha, NE 68114 (402)354-4541 FAX: (402)354-4535

April, 2021

Globulin reference range change

Effective May 4, 2021, Methodist Health System laboratories will change the reference range for globulin.

Background:

After a thorough review of globulin results over the last year in patients in the Methodist Health System it was determined that the reference range for males and females > 12 years of age should be updated.

The previous reference range was: 2.1 to 3.7 g/dL. **The NEW reference range is: 1.7 to 3.7 g/dL**.

Globulin is a generic term for most blood proteins other than albumin. It is a *calculated* value and represents the difference between total serum protein and serum albumin levels.

Order:

Globulin level is part of Comprehensive Metabolic panel (Comp Panel/CMP) and the Hepatic Function Panel (Liver Panel)

Specimen Requirements:

5 mL in SST (serum separating tube with gel separator); PST (heparinized plasma separating tube with gel separator) is also acceptable.

Centrifuge and send SST or serum. Refrigerate 2-8°C.

Refrigerated - 4 days stability. Freeze for extended storage after separation from cells and gel.

Test Performed:

Performed 24/7 with an expected turnaround time (TAT) of 3 hours.

Please direct any questions to Dr. George Bedrnicek at 402 955-5528 or to Dr. Deborah Perry at 402 354-4559.